# Sterile Autoclavable Container for Nisco Encapsulation unit VAR J30

The sterile container for the Nisco J30 is designed for encapsulation with the aerodynamic nozzle. It is the ideal tool if you have to produce beads on high level pharmaceutical or medical standard. It is delivered with the required material certificates to ease the validation of your process.



### Legend:

- 1 Product feed 5
- Air inlet from pressure 6 container with sterile 0.2μm filter
- 3 Vent to hold the vessel in pressureless condition with sterile filter o.2µm
  4 Clamp cover
- nozzle Autoclavable glass vessel (process is visible from outside) Hardening solution

Sterile aerodynamic

- with beads
- Agitation element for stirring

### Working Principle: Aerodynamic Jetting

7

8





## **Nisco Engineering AG**

Wehntalerstrasse 562 CH-8046 Zurich, Switzerland Tel: +41 44 380 06 30 Fax: +41 44 380 06 31 e-mail: mailbox@nisco.ch http://www.nisco.ch

The product enters through a central needle. The exit orifice, which is centrally in line with the axis of the needle, has been counter sunk externally. The counter sunk leads to the aero dynamical effect that the jet has the smaller diameter when it is passing the orifice than the needle. The needle is enclosed in a pressure chamber with an exit through the orifice. The size of the drops is determined by the product flow rate and the pressure inside the chamber.



### Materials in contact with product

All materials are supplied with material certificates.

Stainless steel parts:	1.4435 or equivalent
Electrode:	titanium
Isolation parts:	PEEK
Nozzles:	1.4435 or equivalent
Gaskets:	EPDM or Silicone depending
	on application
Glass cylinder:	borosilicate
O-Rings:	EPDM or Silicone depending
	on application
Hoses:	silicone or polypropylene

### Materials without contact to product

Stainless steel parts:

1.4301 or equivalent