

# Encapsulation Systems VAR J1

## Complete Installation




### Nisco Engineering AG

Dufourstrasse 110  
CH-8008 Zurich, Switzerland

Tel: +41 44 380 06 30  
Fax: +41 44 380 06 31  
e-mail: [mailbox@nisco.ch](mailto:mailbox@nisco.ch)  
<http://www.nisco.ch>

### Legend



1. **Pressure Reducing Station** consisting of precision pressure reducing Valve, manometer and fixation bar for mounting on a support

2. **Suspension Flow Meter** made of aluminium consisting of floating body 0-8 L/min and fixation bar for mounting on a support

3. **Stainless Steel Lab Support** (base plate with bar 12mm diameter with fixation devices for the flow meter, the pressure reducing station and the VAR J1 encapsulation unit)

4. **Encapsulation Unit VAR J1** including an adapter, a sterile filter for air and one set of spare parts (gaskets and spare needle)

5. **Beaker**

6. **Agitator**

7. **Support** for Syringe Pump height adjustment made of aluminium

8. **Syringe Pump** KDS Model 100 (KDS Scientific)